en Instructions for Use PPE Products for Biogel Tech, Biogel NeoTech, Biogel PI Tech, Biogel PI Tech Indicator Underglove.

These gloves are in conformity with Regulation (EU) 2016/425 (Personal Protective Equipment). The gloves comply with EN420:2003+A1:2009.

USE:
For Single Use Only
Regulation (EU) 2016/425 category III product of complex design. The gloves are designed to provide a barrier to microorganisms. The gloves are liquid resistant and can be used for splash protection against certain chemicals. Please consult Mölnlycke Health Care if more information is required on the level of performance of the gloves with specific chemicals.
Chemical data is available from Mölnlycke Health Care.

The actual duration of protection provided in the workplace may vary considerably from these performance levels due to other factors influencing the performance, such as temperature, abrasion, degradation.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.

The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture.

It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.

When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

Before usage, inspect the gloves for any defect or imperfections.


The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PRECAUTIONS:
Gloves damaged in use should be removed and replaced as soon as damage is detected.

DISPOSAL:
Used gloves/sleeves may be contaminated with infectious or other hazardous materials. Dispose of according to Local Authority Regulations.

EU type examination certificate from BSI Group [No 2797] and is the notified body for Regulation (EU) 2016/425: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands

For Declaration of Conformity, please refer to www.molnlycke.com/Declaration-of-Conformities

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
Tested Chemical : 40% Sodium Hydroxide				
No.	Product Description	Product Code	Permeation Performance Level (EN16523-1:2015)	Chemical degradation, Mean % (EN374-4)
1	Biogel PI Tech	44409	Level 6	-58
2	Biogel PI Tech Indicator Underglove	44416	Level 6	-95
3	Biogel NeoTech	44509	Level 6	-19
4	Biogel Tech	44822	Level 6	-66
5	Biogel Tech	44961	Level 6	-108

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Ref no: N/A
Component Type: IFU
Created by: Anders Hansen
Artwork version: 1
Date: 191212

Colour information

 Pantone Green

Mölnlycke Health Care, Gamlestadvägen 3C, SE-402 52, Göteborg.
This artwork was created based on internal procedure T-085.

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Mölnlycke Health Care AB