

# Mölnlycke®


## Avance® Solo Adapt

NEGATIVE PRESSURE WOUND THERAPY

Avance® Solo Adapt Foam, Film and Transfer Port  
Avance® Solo Adapt Film

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
[www.molnlycke.com](http://www.molnlycke.com)

**R<sub>x</sub>** only Caution: Federal [US] Law restricts this device to sale by or on the order of a physician [or properly licensed practitioner].

**MD**








- en** Medical Device
- de** Medizinprodukt
- fr** Dispositif médical
- es** Productos sanitarios
- nl** Medisch hulpmiddel
- sv** Medicinteknisk produkt
- it** Dispositivo medico
- fi** Lääkinnällinen laite
- pt** Dispositivo médico
- no** Medisinsk utstyr
- ar** جهاز طبي

- en** Do not use if package is damaged
- de** Bei beschädigter Verp. verwerfen
- fr** Ne pas utiliser si l'emballage est endommagé
- es** No utilizar si el envase está dañado
- nl** Niet gebruiken indien verpakking beschadigd is
- sv** Skall inte användas om förpackningen är skadad
- it** Non utilizzare se la confezione è danneggiata
- fi** Älä käytä jos pakkaus vaurioitunut
- pt** Não usar se embalagem danificada
- no** Må ikke brukes hvis forpakningen er skadet
- ar** لا تستعمل إذا كان غلاف التغليف الداخلي تالفاً.

# Avance<sup>®</sup> Solo Adapt

NEGATIVE PRESSURE WOUND THERAPY

Patent: [www.molnlycke.com/virtual-patent-marking/](http://www.molnlycke.com/virtual-patent-marking/)

Australian sponsor address:  
Mölnlycke Health Care Pty. Ltd.  
12 Narabang Way, Belrose, NSW 2085 AUSTRALIA

Toll free number:  
USA 1-800-882-4582  
Canada 1-800-494-5134

# Avance® Solo Adapt

## Avance® Solo Adapt Foam, Film and Transfer Port

### Avance® Solo Adapt Film

Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film are products within the Avance Solo Adapt Negative Pressure Wound Therapy (NPWT) System. Avance Solo Adapt NPWT System consists of Avance Solo Adapt Pump, Avance Solo Canister 50 ml, Avance Solo Adapt Foam, Avance Solo Adapt Film and Avance Solo Adapt Transfer Port which together form a system for wound management via the application of negative pressure.

These instructions for use contain information relevant to the Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film for Health Care Professionals. Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film are intended for use by healthcare professionals for therapy on patients in healthcare facilities and home care settings.

For complete list of warnings, precautions and instructions regarding the Avance Solo Adapt NPWT System, including those relevant to the Avance Solo Adapt Pump and Avance Solo Canister see HEALTH CARE PROFESSIONAL Instructions for use provided with Avance Solo Adapt NPWT System.

#### 1. Product description

Avance Solo Adapt Foam, Film and Transfer Port consists of

**Avance Solo Adapt Foam**, a single use polyurethane foam wound filler, distributing negative pressure to the wound bed and facilitating transport of exudate

**Avance Solo Adapt Film**, a single use breathable soft silicone film, to fixate the foam and a to achieve a seal

**Avance Solo Adapt Transfer Port**, a single use polyolefin/polyurethane transfer port, for transporting exudate

Avance Solo Adapt Film consists of

**Avance Solo Adapt Film**, a single use breathable soft silicone film, to fixate the foam and a to achieve a seal

Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film are applicable for use on low to moderately exuding wounds up to the size of the Avance Solo Adapt Foam.

Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film are intended for adults.

#### Material content

*Film:* polyethylene, polyurethane, soft silicone, polyacrylate adhesive

*Foam:* polyurethane

*Transfer Port and tubing with clamp:* Polyurethane, polyester, polyacrylate adhesive, polyolefin based thermoplastic elastomer, polyethylene

*Connectors:* acrylonitrile butadiene styrene copolymer, thermoplastic olefin, polyethylene

#### 2. Indications for use

The Avance Solo Adapt NPWT System is indicated for patients who would benefit from wound management via the application of negative pressure, particularly as the device may promote wound healing through the removal of exudate, infectious material.

Avance Solo Adapt NPWT System is indicated for removal of low to moderate amounts of exudate from pressure ulcers.

#### 3. Contraindications

Avance Solo Adapt NPWT System is contraindicated for patients with the following conditions:

- Malignancy in the wound or margins of the wound
- Untreated and previously confirmed osteomyelitis
- Non-enteric and unexplored fistulas
- Necrotic tissue with eschar present
- Exposed nerves, arteries, veins or organs
- Exposed anastomotic site

#### 4. Warnings

- Application of NPWT may increase the risk of bleeding. If sudden or increased bleeding is observed, immediately take appropriate action to stop bleeding and seek emergency medical attention.
- Patients at high risk of bleeding, such as patients receiving anticoagulant therapy or with altered hemostasis, must be monitored closely during therapy. Ensure to establish hemostasis before applying therapy.
- Patients at risk of bleeding complications due to e.g. a history of vascular anastomosis or friable, irradiated, sutured or infected blood vessels, shall be monitored carefully during therapy.
- Patients undergoing NPWT need frequent supervision. Routinely check that the negative pressure therapy is active, the dressing should be contracted and firm to the touch. If it is necessary to discontinue therapy, the time elapsed without negative pressure must comply with instructions from healthcare professional.
- During therapy, make sure that the pump and tubing are positioned to eliminate the risk of
  - imprints
  - contamination
  - entrapment or strangulation
  - kinking or blocking of tubing
  - being exposed to sources of heat
- Cover or remove sharp edges or bones in the wound prior to the application of the dressing due to the risk of puncturing of organs and blood vessels.
- Tendons, ligaments and nerves shall be protected from direct contact with Avance Solo Adapt Foam.
- Do not place Avance Solo Adapt Foam into unexplored or blind tunnels.
- If defibrillation is required, disconnect Avance Solo Adapt Pump and remove the dressing if the dressing positioning interferes.
- Avance Solo Adapt Pump is Magnetic Resonance (MR) unsafe, do not take into a Magnetic Resonance Imaging (MRI) environment. Avance Solo Adapt Foam, Avance Solo Adapt Film and Avance Solo Adapt Transfer Port are MR safe. The impact of the foam, film and transfer port on Magnetic Resonance Tomography (MRT)/Magnetic Resonance Imaging (MRI) imaging artifacts is unknown.
- If the patient experiences sudden elevation in blood pressure or heart rate in response to stimulation of the sympathetic nervous system, immediately stop the therapy to help minimize sensory stimulation and seek emergency medical attention.
- To minimize the risk of bradycardia, do not apply the dressing close to the vagus nerve.
- Do not use oxidizing agents such as hypochlorite solutions or hydrogen peroxide prior to use of the Avance Solo Adapt Foam.

- Products in Avance Solo Adapt NPWT System contain small parts which could be a potential choking hazard. Keep this device out of the reach of children.
- Keep this device out of the reach of pets.

## 5. Precautions

- Before starting therapy, evaluate the patient's nutritional status and address severe malnutrition. Also consider the patient's size and weight.
- Signs of possible infection or complications must be addressed immediately. Monitor the device, wound, surrounding skin, and patient status accordingly to ensure effective and safe treatment and patient comfort.
- For patients with ischemic condition or at application of a circumferential dressing, extra monitoring of wound status is required to avoid risk of compromised circulation.
- Multiple layers of film or addition of an occlusive dressing over the Avance Solo Adapt Film may decrease the breathability, which may increase the risk of maceration.
- Application of certain skin protection products or using cleansing products prior to the application of Avance Solo Adapt Film, can affect the ability of the film to adhere securely and create sufficient sealing.
- Do not stretch Avance Solo Adapt Film since it may cause damage to the surrounding skin when negative pressure is applied.
- Do not use the products on patient and/or user with known hypersensitivity to the materials/components of the products.
- Ingrowth of tissue may occur if Avance Solo Adapt Foam is not changed in accordance with recommendations or as appropriate for the wound condition of the individual patient.
- Consider use of non-adherent wound contact layer to protect fragile tissue.
- Do not place Avance Solo Adapt Foam onto intact skin or non-protected superficial or retention sutures.
- For daily hygiene routines do not expose the dressing to extensive contact with water.
- Avance Solo Adapt Foam and Avance Solo Adapt Film shall only be applied and changed by healthcare professional.
- No modification of this device (foam, film and transfer port) outside the instructions for use is allowed as modifications may significantly compromise the ability of the system to deliver therapy.
- Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film, shall only be used within the Avance Solo Adapt NPWT System.
- Do not cut the tubing or detach tubing from the canister.
- Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film are provided sterile. Do not use if the sterile barrier is damaged or has been opened prior to use. Do not re-sterilize.
- Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film are single use.
- Do not reuse Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film. If reused, performance of the product may deteriorate and/or cross contamination may occur.
- Always clamp the transfer port tubing and the canister tubing before disconnecting them.

## 6. Instructions for use

### 6.1. To consider before use

Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film are applicable for use on low to moderately exuding wounds up to the size of the Avance Solo Adapt Foam. The maximum allowed wound size area shall not exceed the size of Avance Solo Adapt Foam [67 cm<sup>2</sup>/ 10.4 inch<sup>2</sup>]. Do not use more foam than what is supplied in one Avance Solo Adapt Foam, Film, Transfer port kit.

Low exuding wounds are considered to be up to 0.6 g/ cm<sup>2</sup> wound area/ 24 hours [0.6 g/ 0.16 inch<sup>2</sup> wound area/ 24 hours]. Moderate exuding wounds are considered to be up to 1.1 g/ cm<sup>2</sup> wound area/ 24 hours [1.1 g/ 0.16 inch<sup>2</sup> wound area/ 24 hours]. 1 g exudate is considered equal to 1 ml.

The film shall cover the foam and adhere to the surrounding skin with a margin of 4-5 cm/ 1.6- 2.0 inches. It is preferable to apply the film as a single piece. In case several pieces are required, ensure overlap of the pieces of 4-5 cm/ 1.6- 2.0 inches for proper adhesion. Minimize number of overlaps to avoid moisture build up under the film and/or the risk of air leakage.

If there is a clinical need to treat an underlying pathology of venous insufficiency with compression therapy, ensure that Avance Solo Adapt NPWT System do not induce any pressure points that may result in discomfort or tissue damage to the patient. Do not place the tubing and transfer port under any form of compression garment or bandage.

### 6.2. Pump Set Up

For pump set up see HEALTH CARE PROFESSIONAL Instructions for use provided with Avance Solo Adapt NPWT System.

### 6.3. Dressing Application Procedure

To apply, use clean/aseptic or sterile techniques in accordance with local protocol.

1. Cleanse and debride the wound bed in accordance with local protocol.
2. Cleanse the peri wound skin and pat dry. Apply skin protector/ skin sealant on peri wound skin if required. Certain skin protection products or skin lotions can affect the ability of the film to adhere securely.
3. Eliminate or cover sharp edges and bone fragments with non-adherent wound contact layer due to risk of puncturing organs or blood vessels while under negative pressure. Document the use of wound contact layer in the patient's notes.
4. Cut the foam into an appropriate size corresponding to the dimensions of the wound cavity.
5. Do not cut the foam above the wound site as fragments may fall into the wound. Make sure that no foam fragments are left in the wound or on wound edges.
6. Sufficiently fill the wound cavity without overpacking, as this could damage tissue, affect exudate removal, or affect delivery of negative pressure. Ensure contact between all pieces of foam placed in the wound. Consider use of a non-adherent wound contact layer to protect fragile tissue.
7. Do not place Avance Solo Adapt Foam onto intact skin or non-protected superficial or retention sutures.
8. Document the number of pieces of foam used in the patient's notes.
9. Cut the film to the appropriate size allowing a margin of 4-5 cm/ 1.6-2.0 inches onto surrounding skin.
10. Grip the center part of the release film and pull to expose the adhesive surface.
11. Without stretching, position the film centrally over the wound.
12. Gently remove the remaining release films, smoothing down the film to minimize wrinkles and creases. Do not stretch the film during application as this may cause damage to the surrounding skin when negative pressure is applied.

13. Ensure that the film is in contact with the foam.

14. Contamination of silicone adhesive on top of the film prior to application of the Avance Solo Adapt Transfer Port can affect the ability of the Avance Solo Adapt Transfer Port to adhere securely.

15. Cut a hole, < 2 cm/ 0.8 inch in diameter, in the film covering the wound cavity.

16. Remove the release film on the transfer port and apply it directly over the cut hole in the film. Press gently to fixate it on top of the film. For wounds smaller than the transfer port, protect peri wound skin from imprints and tissue damage, with Avance Solo Adapt Film and Avance Solo Adapt Foam. The tubing should be laid to prevent bends and kinks of the tubing and to avoid risk of imprints and discomfort to the patient.

17. Connect the canister tubing to the transfer port tubing using the connectors on the end of each tubing. Ensure that the tubing is not clamped.

18. To start therapy, press and hold the green pump start button, release after two (2) seconds so that the pump activates and the green start button flashes. Negative pressure shall be achieved within two (2) minutes after pump start.

19. Once the negative pressure is applied, the dressing will contract and be firm to the touch. A wrinkled appearance of the dressing indicates that negative pressure has been achieved and maintained.

20. If negative pressure is hard to achieve, adjust film or press around the edges to improve contact with the skin. Take care to minimize wrinkles and gaps. Consider using additional Avance Solo Adapt Film or start over and reseal, if necessary.

### 6.4. Frequency of dressing change

Avance Solo Adapt Foam, Film and Transfer Port should be changed every 48 to 72 hours, but no less than 3 times a week, or as instructed by healthcare professional.

NOTE: Ingrowth of tissue may occur if the foam is not changed in accordance with recommendations or as appropriate for the wound condition of the individual patient.

### 6.5. Dressing removal procedure

To remove the dressing, perform the following steps:

21. If the pump is active, pause therapy by pressing down the green start button, release after two (2) seconds.
22. Block both the canister tubing and transfer port tubing by positioning the slide clamps next to the connector and slide across the tubing until secured.
23. Disconnect the canister tubing from the transfer port tubing by squeezing the connector from both sides and pull apart.
24. Gently peel back one corner of the film and stretch to facilitate breakage of the seal. Proceed with this technique (in the direction of hair growth) until the film is completely removed.
25. Remove the foam gently. If adherence of foam to wound bed is observed, consider moistening the foam. Make sure no fragments are left in the wound.

If the patient experiences pain during dressing removal, consider the use of pain relief medication, as instructed by a healthcare professional, when changing dressings. Consult patient's notes to ensure that all materials used have been removed.

## 6.6. Changing the canister, changing the batteries and daily use

For instructions how to change the canister or the batteries and for daily use see HEALTH CARE PROFESSIONAL Instructions for use provided with Avance Solo Adapt NPWT System.

## 6.7. Disposal

Avance Solo Adapt Foam, Film and Transfer Port and Avance Solo Adapt Film are single use products. After use, dispose of the products as clinical waste in accordance with local regulations.

For more information about safe disposal, see [www.mölnlycke.com/wastehandling](http://www.mölnlycke.com/wastehandling) or contact your local Mölnlycke Health Care representative.

## 7. Avance Solo Adapt NPWT System Indicators, Alarms and Troubleshooting

For system indicators for normal use, alarms and troubleshooting HEALTH CARE PROFESSIONAL Instructions for use provided with Avance Solo Adapt NPWT System.

## 8. Caution

Avance Solo Adapt Foam, Film and Transfer Port and single packed Avance Solo Adapt Film must be used in accordance with these Instructions for Use in combination with the HEALTH CARE PROFESSIONAL Instructions for use provided with Avance Solo Adapt NPWT System. Read the instructions before using the system and have them available during use. Failure to read and understand these instructions may lead to misuse of the system and improper performance. These instructions are a general guide for the use of the product. Specific medical situations must be addressed by a healthcare professional.

## 9. Other information

If any serious incident has occurred in relation to Avance Solo Adapt Foam, Film and Transfer Port or Avance Solo Adapt Film, it should be reported to Mölnlycke Health Care.

Avance® and Safetac® are registered trademarks of Mölnlycke Health Care AB

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